



Do urinary BMP levels change with disease progression in idiopathic and genetic PD?

Kalpana Merchant, TransThera Consulting Co.
Frank Hsieh, Nextcea, Inc.

PPMI Project ID: 145

Summary

Urine biospecimens from the LRRK2 genetic cohort of PPMI have been analyzed for three phospholipids: total di-BMP-22:6, its distinct 2, 2'-isoform (2, 2' di-22:6-BMP) and total di-BMP-18:1. The following table lists the biospecimens analyzed.

Study Arm	BL Visit	V04 (1 year)	V06 (2 year)	
LRRK2 PD	150	121	82	
LRRK2 Unaffected	202	159	100	
Total	352	280	182	866

Method

Analysis of BMP isoforms:

Quantitative Ultra Performance Liquid Chromatography Tandem Mass Spectrometry (UPLC-MS/MS) analyses: Urinary BMPs were extracted by liquid-liquid extraction and measured by targeted UPLC-MS/MS as described previously by Liu et al., 2014. Absolute levels of total di-22:6-BMP (the sum of its three isoforms), its 2, 2'-isoform (2,2'-di-22:6-BMP), and total di-18:1-BMP were measured using a SCIEX X500B QTOF LC-MS/MS system (SCIEX, Framingham, MA). Injections were made using a Shimadzu Nexera XR UPLC system (Shimadzu Scientific Instruments, Japan). The instruments were controlled by SCIEX OS Software. Quantitation was performed using authentic di-22:6-BMP and di-18:1-BMP reference standards. Di-14:0-BMP was used as an internal standard.

Calibration and data processing: The intensities of the analytes and internal standards were determined by integration of extracted ion peak areas using SCIEX OS-Q Software. Calibration curves were prepared by plotting the peak area ratios for each analyte to internal standard versus concentration. The model for the calibration curve was linear with $(1/x^2)$ weighting. Measured concentrations of the urinary lipids (ng/mL) were normalized by the concentration of urine creatinine and reported in ng/mg creatinine.





Creatinine measurement:

Concentrations of urinary creatinine were measured by colorimetric assay (method of Jaffé) with Parameter Creatinine Assay test reagents (R&D Systems, Minneapolis, MN) using a BioTek ELx800 absorbance microplate reader with Gen5 Microplate Reader and Imager Software 2.09 (Fisher Scientific, Hampton, NH).

References

Liu N, Tengstrand EA, Chourb L, Hsieh FY. Di-22:6 bis(monoacylglycerol)phosphate: A clinical biomarker of drug-induced phospholipidosis for drug development and safety assessment. Toxicology and applied pharmacology 2014;279:467-476).

About the Authors

This document was prepared by Kalpana Merchant (TransTHera Consulting Co.) and Frank Hsieh (Nextcea, Inc.). For more information please contact Kalpana Merchant at merchant.transthera@gmail.com or Frank Hsieh at frank.hsieh@nextcea.com

Notice: This document is presented by the author(s) as a service to PPMI data users. However, users should be aware that no formal review process has vetted this document and that PPMI cannot guarantee the accuracy or utility of this document.

